


## LABS Monitoring & MINIZ Program



## **Features**

- Check the temperature, liquid nitrogen level, alarm conditions
- Screen Display Temperature, Liquid Nitrogen Level, Alarm conditions, Liquid Nitrogen supply condition, and equipment number & quantity of the alarmed units
- Saved data
- · Controls supply of the liquid nitrogen
- One fill All fill operation
- Wireless WIFI Comunication
- · All of the details can be easily seen at the FULL HD monitor
- Connected with O2 monitor system that it shows the Oxygen density of the room where cryopreservation equipments are installed so it safeguards whether worker is safe to enter the room before going in.

## **Specification**

- LAMINIZ program
- Program Operating System : WINDOW XP or upper version
- Communication Method : Wireless WIFI Comunication
- No. of equipments possible for connection : Maximum 150 units

25

Computer Specification:

RAM: DDR3 / 8GB or above

CPU: Intel i5 or above

SSD: 128GB or above

ODD: DVD-RW or above

OS: Windows 10 pro 이상

Connection type : Wireless way (LAN Port 연결)

Converter: Wifi Data converter

Monitor: Full HD Class

- Solenoid Valve Controller
- · Remote Alarm Contact Controller